Hersteller / Manufacturer

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
GB

Intended Use

The Papillon® nasal mask is designed for the treatment of patients suffering from restrictive, obstructive, or central respiratory disorders who require nasal positive airway pressure therapy (nCPAP or biLEVEL). Examples of these conditions include sleep apnea, chronic restrictive or obstructive lung diseases, and deformations of the rib cage.

restrictive = limiting
obstructive = blocking
central = brain-triggered

The nasal mask may be used by one patient only, and may not be made available for use by another patient.

The Papillon nasal mask is a high-quality MAP Medizin-Technologie GmbH product, and has been developed on the basis of leading-edge research into respiratory therapies.

Legend

See Figures on inside front cover

- 1 Air tubing
- 2 Tubing connector
- 3 Angled adapter
- 4 Connectors for external pressure measurement and for administering O₂
- 5 Connector cover
- 6 Mounting holes for the forehead pad
- 7 Forehead frame
- 8 Forehead pad
- 9 Slots to attach the headgear
- 10 Mask frame
- 11 Frame base
- 12 Mask cushion
- 13 Exhalation vents
- 14 Tongue joint
- 15 Groove joint
- 16 Scale showing the four hinge positions
- 17 Hinge grip
- 18 Clip fastener
- 19 AeroFix® plus tubing clip (not included, REF 632 006)
- 20 AeroFix® II headgear

Before Using for the First Time

Please read and note carefully all the information and instructions provided in the Instruction Manual. Hygienic use of the Papillon nasal mask and product durability in accordance with the useful life specifications can be achieved only by strict compliance with these instructions. Do not use the Papillon nasal mask if the packaging is damaged. Check that the complete package contents are present. The package contents comprise the nasal mask, the AeroFix® II headgear [20], the clip fastener [18], and this Instruction Manual. Clean all components and the headgear of the new mask before using them for the first time. See the "Disassembly of the Nasal Mask" and "Cleaning and Care" sections for more information. Wash the headgear separately for the first few times, as the color may run.

Safety Information

- The type of mask used affects the therapeutic pressure. The treating physician should therefore check whether the pressure needs to be adjusted if you were previously being treated with a different mask.
- Please contact your treating physician if you experience pain in the nose, sinuses, or ears after commencing treatment, or if you have an allergic reaction to any of the materials used.
- The nasal mask is suitable for pressures ranging between 4 and 20 hPa.
- To prevent exhaled air being rebreathed, it is essential that this air is always able to escape through the exhalation vents [13] in the mask cushion [12]. Before each use, always ensure that the exhalation vents are not blocked. Any deposits should be removed.
- When fitting spare parts, the colour of the mask cushion is important. Only ice-blue cushions should be used with the Papillon mask. The ice-blue mask cushion with integrated exhalation vents used for the Papillon nasal mask must on no account be interchanged with the mint-green mask cushion of the Silent Papillon nasal mask.
- The connection system for the Papillon nasal mask must never be interchanged with the exhalation system of the Silent Papillon nasal mask.

- Avoid excessive build-up of condensed water, particularly at low therapeutic pressures. If necessary, reduce the settings of the air humidifier accordingly.
- You should remove the mask as soon as you switch off the respiratory therapy device or detach the mask from the air tubing. This will prevent you from rebreathing your own exhaled air.
- To remove the mask safely, first open the clip fastener [18] on the headgear [20] and then detach the tubing connector [2] from the angled adapter [3]. Do not pull on the mask.
- Follow the instructions on regular cleaning of the mask to prevent possible respiratory tract infections.
- The Papillon nasal mask may be used only for nasal positive airway pressure therapy (nCPAP or biLEVEL).
- Do not use or repair a defective mask. Use original spare parts only.

Fitting the Nasal Mask

Always wash your face and hands before using the mask.

- [Fig. C] Insert the tubing connector [2] into the air tubing [1].
- [Fig. D] Thread the strap of the clip fastener [18] through one of the lower slots [9] for the headgear.
- Lay out the AeroFix II headgear [20] in front of you, so that the side with the washing instructions is facing away from you. The longer straps should be at the top and the shorter ones at the bottom.
- [Fig. E] Open the Velcro fasteners on the AeroFix II headgear and thread the ends through the appropriate slots [9] on the mask frame and clip fastener. Always thread from the inside and close the Velcro fasteners on the outer side of the straps.
- [Fig. E] Optional: attach the middle of each upper strap with the AeroFix® plus tubing clip [19], available separately.
- [Fig. F] Stand in front of a mirror and hold the nasal mask on your face using one hand. With the other hand, pull the headgear over your head with the clip fastener open.
- [Fig. G] Close the clip fastener and readjust the ends of the headgear one by one so that they are even. Ensure that the upper straps are at the same height above each ear. The lower straps must be below the ears. The

saddle part of the headgear should be centered at the back of your head.

- [Fig. H] Push the tubing connector together with the air tubing onto the angled adapter [3]. Optional: then secure the air tubing to your head using the Velcro fastener of the AeroFix plus tubing clip.
- Switch your respiratory therapy device on and get into bed, in order to make sure the mask is fitted correctly.
- The mask should not pinch, and should not leave any pressure marks. No air should be able to escape between the skin and the mask cushion. Use the headgear to adjust the fit of your Papillon nasal mask accordingly.
- [Fig. I] It is also possible to adjust the angle between the frame base [11] and the forehead frame [7]. There are four possible positions [16]. To make the adjustment, take off the mask and squeeze the hinge grip [17] with one hand while adjusting the position of the forehead frame with the other hand. If necessary, you can also attach the forehead pad [8] onto the upper mounting holes [6] of the forehead frame.

Disassembly of the Nasal Mask

- [Fig. A] Grip the tubing connector [2] and detach it from the angled adapter [3] by snapping it off to the side.
- Push the mask cushion [12] backwards out of the frame base [11].
- Remove the forehead pad [8] from the forehead frame [7].
- Detach the angled adapter from the mask frame [10] by twisting the adapter anti-clockwise, unlocking it, and pulling it out.
- Remove the connector cover [5] from the angled adapter.
- Detach the forehead frame from the frame base by squeezing the hinge grip [17] while pulling the forehead frame out backwards along the guide rail.

Assembly

- [Fig. A] Squeeze the hinge grip [17] and push the forehead frame [7] into the guide rail of the frame base [11].
- [Fig. A] Push the connector cover [5] onto the connectors for external pressure measurement and administering oxygen [4] on the angled adapter [3].

c) [Fig. B-a] Push the angled adapter from above into the round opening of the frame base. Insert the two different-sized catches into the corresponding slots. Lock the angled adapter by turning it clockwise to the limit stop.

d) [Fig. A] Push the mask cushion [12] from the rear through the frame base until it is below the angled adapter.

Push the round opening on the upper part of the mask cushion as far as possible over the angled adapter. There is a groove below the six exhalation vents [13]. This groove must lock onto the frame [Fig. B-b]. Be careful to ensure that the tongue joints [14] snap into the corresponding groove joints [15] on the sides of the mask cushion [Fig. B-c].

e) [Fig. A] Press the two pegs on the forehead pad [8] from the rear into the lower mounting holes [6] of the forehead frame.

Cleaning and Care in the Home Environment

For reasons of hygiene, the nasal mask must be cleaned every day.

Wash your hands before cleaning the nasal mask, and always clean and rinse the mask with drinking-quality water. Never use soap, chlorinated, corrosive, or abrasive products, or cleaning agents with a high alcohol content, since these may damage the material. Never use metal brushes or steel wool.

Cleaning

Clean all components of the nasal mask daily after use with a warm (approx. 30 °C), mild dishwashing solution for about 5 minutes. For effective cleaning, take the mask apart as described in the “Disassembly of the Nasal Mask” section. Remove any deposits from the exhalation vents [13]. For optimum cleaning, the mask cushion can be turned inside out.

Then rinse all parts thoroughly in drinking-quality water to remove all remaining traces of dishwashing solution.

The AeroFix II headgear [20] and AeroFix plus tubing clip [19] should be cleaned once weekly. Detach the headgear and tubing clip from the nasal mask and close the Velcro fasteners.

The headgear may be washed by hand or in the washing machine at 30 °C. Dry the headgear on slow spin or drip dry.

NB: The headgear is sensitive to heat, and may not be ironed!

Drying

For reasons of hygiene, it is important to dry the nasal mask thoroughly. Wipe the cleaned components dry with a clean, lint-free cloth. Ensure that the nasal mask is not exposed to direct sunlight or placed on a heater, since this could damage the mask material.

Thermal Disinfection

Disinfection of the nasal mask in the home environment is generally not required. Disinfection of the mask may be necessary where infectious diseases are present, in the case of a known immune deficiency, or where directed by a physician. The purpose of the disinfection process is to kill pathogens. For disinfecting your nasal mask, we recommend the use of a commercial baby bottle disinfectant. The nasal mask can be thermally disinfected up to 60 times using such a device.

Warning! Burns hazard. Always comply strictly with the contents of the disinfectant instruction manual.

Before disinfection, the nasal mask must always be properly cleaned, rinsed, dried, and fully dismantled [see under “Disassembly of the nasal mask,” “Cleaning,” and “Drying”].

Always disinfect with distilled water, to prevent any lime deposits from building up on the surface of the nasal mask. The water vapor must reach all parts of the nasal mask. Allow the nasal mask to cool before removing it from the baby bottle disinfectant. After disinfection, the nasal mask must be dried once again [see under “Drying”].

Functional Test

The correct functioning of the nasal mask must be tested after each cleaning or disinfection operation. Inspect the nasal mask for any defects or changes in the condition of the materials. A defective mask must always be replaced.

Assuming the cleaning procedure described above is followed, the Papillon nasal mask, AeroFix II headgear, and AeroFix plus tubing clip have a useful life of one year.

Storage

Keep the nasal mask in a dry, dust-free location protected from direct sunlight until next required.

Cleaning and Disinfection in the Hospital Environment

Information on cleaning and disinfection procedures for the Papillon nasal mask in the hospital environment is available on request from MAP Medizin-Technologie GmbH (REF 560 092).

Disposal

The Papillon nasal mask, AeroFix II headgear, and AeroFix plus tubing clip contain no toxic substances. Dispose of the nasal mask in accordance with applicable national laws and regulations. Except where stipulated otherwise, before being taken to a recycling center the nasal mask should be cleaned in accordance with the procedures described in this Instruction Manual.

Warranty

A two-year warranty applies for manufacturing defects. The warranty does not cover normal wear and tear. For further details, see the MAP Medizin-Technologie GmbH General Terms and Conditions of Business.

MAP Medizin-Technologie GmbH shall not be liable for any impairments to the safety, reliability, or performance of the product caused by:

- use other than that specified in the Instruction Manual
- use in conjunction with other products not designed for this purpose
- any technical alterations made to the product.

The mask must always be cleaned and disinfected before being returned for servicing. On completion of the above procedure, dry the mask and pack it in a securely fastened transport bag.

Symbols

See Instruction Manual

Warning! Failure to follow this instruction could pose a health risk.

REF Part number

Lot number

Date of manufacture (month, year)

Washing temperature 30 °C

Do not use chlorine bleach

Do not iron

Do not spin-dry

Technical Data

Standards

CE labeling in accordance with EC directive 93/42/EEC, class IIa

Quality Management: development, manufacture, sale, and service certified in accordance with ISO 13485:2003 quality management principles.

EN ISO 17510-2

(Sleep apnea respiratory therapy, Part 2: masks and application accessories)

Validated Cleaning and Disinfection

The general suitability of the Papillon® nasal mask for effective cleaning and disinfection by the procedures specified in this Instruction Manual has been established by an independent accredited test laboratory.

Flow Resistance

Pressure drop measured with mask cushion size 2
at 50 l/min = max. 0.1 hPa
at 100 l/min = max. 0.4 hPa

Acoustic Pressure

30.8 dB(A), measured with artificial head at 10 hPa according to EN ISO 17510-2

Tubing Connection

Ø 22 mm for elastomer tubing

Weight

Approx. 92 g

Therapeutic Pressure

4 to 20 hPa

Materials Used


Mask (latex-free): silicone, PC
Headgear (neoprene-free): nylon, lycra, spandex, polyurethane, polyester

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© 2005 Copyright by MAP Medizin-Technologie GmbH
The instructions in this manual may differ from the product. This does not affect the product's performance. This instruction manual is protected by copyright. All rights reserved. Any reproduction including excerpts requires the express consent of MAP Medizin-Technologie GmbH.